

Comprehensive Skin Lesion Analysis Report

Prediction: nv


Confidence: 0.84

Generated on: 2025-03-19 08:13:11

This report is for research purposes only.

Image and Predictions

Original Image



Full Prediction Confidence






Image Name: sample_image
Date Analyzed: 2025-03-19 08:13:12
Model: SkinLesionModel
Image Size: 224x224
Top Predictions:
1. nv: 0.8362
2. mel: 0.1554
3. bkl: 0.0048

Explanation Methods - Part 1


Captum Grad-CAM Explanation



Grad-CAM Heatmap




LIME Explanation




ations) perturbs the input image by segmenting it and the model to approximate how the segments affect the model's prediction. Green regions contribute positively to the prediction, while red regions negatively contribute.

Explanation Methods - Part 2


Integrated Gradients Explanation




IG Attribution Map



Class Activation Mapping



CAM Heatmap



the model's decision-making process, but they are approximations and should be interpreted with caution. **Important Information**

For any concerns about skin lesions, please consult a qualified dermatologist. Early detection and proper medical assessment are crucial for skin cancer diagnosis and treatment.

About the Model:

Architecture: SkinLesionModel

Training Dataset: HAM10000

Classes: Melanocytic nevi (nv), Melanoma (mel), Benign keratosis (bkl), Basal cell carcinoma (bcc), Actinic keratoses (akiec), Vascular lesions (vasc), Dermatofibroma (df)

Explanations generated using multiple XAI (Explainable AI) techniques to provide a comprehensive understanding of the model's prediction.